



A 4-Week, Two-Center, Open-Label, Single-Arm Study to Evaluate the Safety and Efficacy of EOPatch in Well-Controlled Type 1 Diabetes Mellitus

Jiyun Park¹, Nammi Park², Sangjin Han³, You-Bin Lee¹, Gyuri Kim¹, Sang-Man Jin¹, Woo Je Lee⁴, Jae Hyeon Kim^{1,5}

¹Division of Endocrinology and Metabolism, Department of Medicine, Samsung Medical Center, Sungkyunkwan University School of Medicine, Seoul, Divisions of ²Clinical Affairs, ³Hardware Development, EOFLOW Co. Ltd., Seongnam,

⁴Department of Internal Medicine, Asan Medical Center, University of Ulsan College of Medicine, Seoul,

⁵Department of Clinical Research Design and Evaluation, Samsung Advanced Institute for Health Sciences & Technology, Sungkyunkwan University, Seoul, Korea

This study evaluated the safety and efficacy of tubeless patch pump called EOPatch in patients with well-controlled type 1 diabetes mellitus (T1DM). This 4-week, two-center, open-label, single-arm study enrolled 10 adult patients diagnosed with T1DM with glycosylated hemoglobin less than 7.5%. The co-primary end points were patch pump usage time for one attachment and number of serious adverse events related to the patch pump. The secondary end points were total amount of insulin injected per patch and changes in glycemic parameters including continuous glucose monitoring data compared to those at study entry. The median usage time per patch was 84.00 hours (interquartile range, 64.50 to 92.50). Serious adverse events did not occur during the trial. Four weeks later, time in range 70 to 180 mg/dL was significantly improved ($70.71\% \pm 17.14\%$ vs. $82.96\% \pm 9.14\%$, $P=0.01$). The times spent below range (<54 mg/dL) and above range (>180 mg/dL) also improved (All $P<0.05$). Four-week treatment with a tubeless patch pump was safe and led to clinical improvement in glycemic control.

Keywords: Diabetes mellitus, type 1; Glycated hemoglobin A; Insulin; Insulin infusion systems

INTRODUCTION

Multiple daily injection (MDI) therapy and continuous subcutaneous insulin infusion (CSII) with an external pump are standard methods of intensive insulin therapy for diabetes. CSII increases patient satisfaction and is effective for glycemic control, prevention of hypoglycemia, and decreasing cardiovascular mortality compared to MDI [1-4]. An insulin infusion set (IIS) delivers insulin from a pump to subcutaneous tissue and consists of a plastic tube connected to a cannula or needle [5]. Problems related to IIS include kinking, dislocation or leakage of cannula, reservoir leakage, a loose connection between the reservoir and cannula, or occlusion [6]. IIS failure is

important because it is the most common problem in patients using insulin pumps, and defects can cause serious adverse effects such as diabetic ketoacidosis (DKA). Patch pumps generally are small and attach directly to the skin without tubing to reduce IIS-related problems [7]. They also are easy to train patients to use, maintain insulin at similar temperatures inside the patch pump, and show low variation in insulin delivery according to position [8].

A Korean company, EOFLOW (Seongnam, Korea), has developed the EOPatch, a disposable patch type insulin pump that was approved by the National Institute of Food and Drug Safety Evaluation (NIFDS) in Korea in December 2017 (Product license No. 17-959). A new model that was upgraded with

Corresponding author: Jae Hyeon Kim <https://orcid.org/0000-0001-5001-963X>
Division of Endocrinology and Metabolism, Department of Medicine, Samsung Medical Center, Sungkyunkwan University School of Medicine, 81 Irwon-ro, Gangnam-gu, Seoul 06351, Korea
E-mail: jaehyeon@skku.edu

Received: Oct. 27, 2021; Accepted: Jan. 5, 2022

This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (<https://creativecommons.org/licenses/by-nc/4.0/>) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

an interlinking blood glucose meter, alarm, and usage time was developed and approved by the NIFDS in June 2020.

We investigate patch pump usage time for one attachment, serious adverse events, and glycemic outcomes in patients with well-controlled T1DM.

METHODS

Study design and procedures

We performed this study at Samsung Medical Center (SMC) and Asan Medical Center (AMC), Korea, from March 2020 through July 2020. The study flow is outlined in Supplementary Fig. 1, and specific training programs are described in the Supplementary Material. Participants used Dexcom G5 to check glucose and advanced diabetes manager (ADM) to control insulin dose and infusion rate (Supplementary Figs. 2 and 3). Bluetooth-connected EOBridge (Supplementary Fig. 4), which is an application for Android and iOS, allows participants and guardians to monitor insulin dose and patch pump operation. In the 4th week, investigators collected the ADM and assessed exact insulin doses and usage times. We also conducted patch pump satisfaction survey (Appendix 1).

Patients

We enrolled adults (19 years of age or older) who were diagnosed with T1DM more than 6 months prior, had used MDI or an insulin pump for more than 3 months, and were in a well-controlled diabetic state (glycosylated hemoglobin [HbA1c] less than 7.5%). A full list of the inclusion and exclusion criteria is provided in the Supplementary Material. Among 10 participants, one dropped out within 3 days due to personal circumstances. The nine remaining patients were included in final analysis. The subjects voluntarily participated in the study and provided written informed consent. The study protocol was conducted in accordance with the principles of the Declaration of Helsinki and Good Clinical Practice guidelines. The protocol was approved by the Institutional Review Boards of SMC (no. 2019-11-220-011) and AMC (no. S2019-2494-0001).

End points

The co-primary end points were patch pump usage time for one attachment and number of serious adverse events. Serious adverse events were defined as severe hyperglycemia (glucose >400 mg/dL for at least 20 minutes), severe hypoglycemia

characterized by altered mental and/or physical functioning that requires assistance from another person for recovery, or DKA [9,10]. The secondary end points were total amount of insulin injected per patch and changes in HbA1c, glycoalbumin, and continuous glucose monitoring (CGM) data.

Statistical analysis

Data are presented as mean \pm standard deviation or median (interquartile range [IQR]) for continuous variables or as number with percentage for categorical variables. Comparisons of glycemic outcomes and CGM data before and after using the investigational device were performed using the paired *t*-test or Wilcoxon's signed rank test. All *P* values were two-tailed, and *P* values <0.05 were considered statistically significant. All statistical analyses were performed using SAS version 9.4 (SAS Institute, Cary, NC, USA).

RESULTS

Baseline characteristics of the study sample are shown in Table 1. The mean age was 38 years, most participants were female

Table 1. Baseline characteristics of the study sample

Patient characteristic	Value
Number	10
Age, yr	38 \pm 12
Male sex	1 (10)
Body mass index, kg/m ²	22.2 \pm 3.5
Years since diagnosis	8.1 \pm 7.7
Duration of insulin use, yr	8.1 \pm 7.7
Conventional pump user	1 (10)
Real-time CGM before trial	
Dexcom	6 (60)
Freestyle Libre	3 (30)
Guardian	1 (10)
Insulin brand	
Humalog (Eli Lilly)	10 (100)
Daily total insulin dose, units	36.90 \pm 11.25
HbA1c, %	6.56 \pm 0.37
Fasting plasma glucose, mg/dL	145.30 \pm 30.78
Current smoker	0

Values are presented as mean \pm standard deviation or number (%). CGM, continuous glucose monitoring; HbA1c, glycosylated hemoglobin.

(90%), and the mean duration of DM was 8.1 years. The mean of total patch usage time is 683.2 hours and the patch was attached for 99.3% of 687.7 hours of study period except for replacement time. The median patch usage time was 84.00 hours (IQR, 64.50 to 92.50). Nine of 84 patches (10.7%) were used for less than 24 hours, six patches (7.1%) were used for less than 48 hours, nine patches (10.7%) were used for less than 72 hours, 11 patches were used for 84 hours, and 49 patches were used over than 84 hours. The reasons for less than 84 hours were due to inlet blockage, pump malfunction, or ADM smart controller malfunction, unexplained hyperglycemia that occurs without apparent medical, dietary insulin dosage, or pump explanation, and last visit. The median insulin amount per patch was 122.90 units (IQR, 86.05 to 149.86), with a minimum of 2.60 units and maximum of 186.60 units. Serious adverse events related to investigational devices did not occur in any participants. There were no significant differences in HbA1c, glucose management indicator (GMI), and glycoalbumin (Table 2). The time in range (TIR, 70 to 180 mg/dL) significantly increased ($70.71\% \pm 17.14\%$ vs. $82.96\% \pm 9.14\%$, $P=0.01$) (Table 2). The times spent below range (TBR, <54 mg/dL) and above range (TAR >180 mg/dL) also improved (Table 2). The total daily dose of insulin (TDD) was not different in the nine participants before and after using a patch pump, 38.33 ± 10.92 and 37.55 ± 8.99 units, respectively ($P=0.639$). There were significant differences in simple satisfaction scores between conventional pump or MDI and investigational device (5.78 ± 1.99 vs. 8.78 ± 1.09 , $P=0.002$). In the EOPatch satisfaction questionnaire, most participants reported they were satisfied with the device and its usefulness (Supplementary Table 1).

DISCUSSION

We observed a median patch usage time of 84.00 hours (IQR, 64.50 to 92.50), and the median insulin amount per patch was 122.90 units (IQR, 86.05 to 149.86). Serious adverse events did not occur during the trial. After the trial, TIR 70 to 180 mg/dL, TBR <54 mg/dL, and TAR >180 mg/dL were significantly improved. In addition, patients were more satisfied with tubeless patch pump than with MDI or other conventional pumps.

The Omnipod, manufactured by the Insulet Corporation, is the most widely used patch pump, can be used for 3 days (72 hours). Treatment outcomes using the Omnipod system include improved glycemic control and reduction of hypoglycemic events compared to MDI usage and were not inferior to

Table 2. Glycemic outcomes before and after EOPatch usage

Variable	Before	After	P value
HbA1c, %	6.56±0.37	6.39±0.26	0.10
Glycoalbumin, %	18.43±2.89	17.57±2.18	0.14
Mean fasting glucose, mg/dL	145.30±30.78	135.50±14.20	0.17
Glucose management indicator, %	6.83±1.06	6.38±0.52	0.19
TIR 70–180 mg/dL, %	70.71±17.14	82.96±9.14	0.01
TBR <70 mg/dL, %	3.17±1.59	2.03±1.71	0.13
TBR <54 mg/dL, %	1.02±0.87	0.27±0.21	0.02
TAR >180 mg/dL, %	18.03±9.69	12.97±7.92	0.02
TAR >250 mg/dL, %	7.03±10.14	1.78±1.54	0.07

Values are presented as mean ± standard deviation.

HbA1c, glycosylated hemoglobin; TIR, time in range; TBR, time below range; TAR, time above range.

those of conventional CSII [7,11]. The EOPatch uses Teflon cannulas like the Omnipod, minimizing the inconvenience caused by needles [12,13]. The motor mechanism is different from that of actuation of the Omnipod, which is a shape memory alloy [8]. The EOPatch uses an electroosmotic pump that allows precise automatic insulin delivery using electrochemistry, which allows a more compact size and enables correct insulin delivery capacity with less noise and heat than existing pumps [14]. EOPatch usage time is 3.5 days compared to 3 days for the Omnipod. In this study, 71.4% of patch pumps were used for more than 72 hours, and median patch usage time was 84 hours. Based on that, the cost burden should be decreased and allow patients to replace the pump only twice per week with regular usage.

This was a short 4-week study that did not allow a sufficient length of time to evaluate HbA1c. Therefore, we observed no significant statistical differences in HbA1c, GMI, or mean plasma glucose level. However, results of CGM data were improved. In particular, mean TIR significantly improved from 70% to 83% without change in TDD, which is a considerable advantage considering that 70% and 80% of each TIR corresponds to HbA1c 7.0% and 6.5% [15]. Considering that all the participants in our study already had used CGM before study entry and HbA1c was less than 7.5%, effective carbohydrate counting and precise insulin dose calculation using bolus calculator might allow glycemic improvements. However, these advantages are not limited to tubeless patch pump and can be experienced with all insulin pumps.

According to T1DM exchange clinical registry data, compared with the fair/poor control group (HbA1c same or over 8.5%), the excellent control group (HbA1c less than 6.5%) was characterized by fewer missing insulin doses and more frequent meal or correction bolus [16]. In other words, with more frequent bolus insulin injections, patients achieve better control of diabetes, highlighting the importance of pump usage. However, CSII is used in only about 5% in T1DM in Korea compared to 63% in United States [17], which can be explained by lack of reimbursement for prescription, education and training of insulin pump by National Health Insurance, and the low prevalence of T1DM [18].

This study has several limitations. First, this is a single-arm study, with a small number of participants and a short observation period. Second, the effects of selection bias could not be excluded because of the nonrandomized study design. Third, there might be volunteer bias because we enrolled patients with well-controlled T1DM (HbA1c less than 7.5%) and provided structured education. Therefore, it might be difficult to generalize our glycemic outcomes to real clinical settings.

In conclusion, the median usage time of the tubeless patch pump was 3.5 days, and 4-week usage was safe and led to clinically meaningful improvement in glycemic control, including severe hypoglycemia and TIR in patients with well-controlled T1DM.

SUPPLEMENTARY MATERIALS

Supplementary materials related to this article can be found online at <https://doi.org/10.4093/dmj.2021.0299>.

CONFLICTS OF INTEREST

Nammi Park and Sangjin Han are employees of EOFLOW. No other potential conflicts of interest relevant to this article were reported.

AUTHOR CONTRIBUTIONS

Conception or design: N.P., W.J.L., J.H.K.

Acquisition, analysis, or interpretation of data: J.P., N.P., S.H., Y.B.L., G.K., S.M.J.

Drafting the work or revising: J.P., N.P., J.H.K.

Final approval of the manuscript: J.P., N.P., J.H.K.

ORCID

Jiyun Park <https://orcid.org/0000-0002-2402-1979>

Jae Hyeon Kim <https://orcid.org/0000-0001-5001-963X>

FUNDING

Funding and material supports for this study were provided by EOFLOW.

ACKNOWLEDGMENTS

We appreciated the trial participants for taking part in the study.

REFERENCES

1. Steineck I, Cederholm J, Eliasson B, Rawshani A, Eeg-Olofsson K, Svensson AM, et al. Insulin pump therapy, multiple daily injections, and cardiovascular mortality in 18,168 people with type 1 diabetes: observational study. *BMJ* 2015;350:h3234.
2. Beato-Vibora P, Yeoh E, Rogers H, Hopkins D, Amiel SA, Choudhary P. Sustained benefit of continuous subcutaneous insulin infusion on glycaemic control and hypoglycaemia in adults with type 1 diabetes. *Diabet Med* 2015;32:1453-9.
3. Hirsch IB, Bode BW, Garg S, Lane WS, Sussman A, Hu P, et al. Continuous subcutaneous insulin infusion (CSII) of insulin aspart versus multiple daily injection of insulin aspart/insulin glargine in type 1 diabetic patients previously treated with CSII. *Diabetes Care* 2005;28:533-8.
4. Pala L, Dicembrini I, Mannucci E. Continuous subcutaneous insulin infusion vs modern multiple injection regimens in type 1 diabetes: an updated meta-analysis of randomized clinical trials. *Acta Diabetol* 2019;56:973-80.
5. Zhang JY, Shang T, Chattaraj S, Cohen O, Heschel M, Vigersky RA, et al. Advances in insulin pump infusion sets symposium report. *J Diabetes Sci Technol* 2021;15:705-9.
6. Deiss D, Adolfsson P, Alkemade-van Zomeren M, Bolli GB, Charpentier G, Cobelli C, et al. Insulin infusion set use: European perspectives and recommendations. *Diabetes Technol Ther* 2016;18:517-24.
7. Lebenthal Y, Lazar L, Benzaquen H, Shalitin S, Phillip M. Patient perceptions of using the OmniPod system compared with conventional insulin pumps in young adults with type 1 diabetes. *Diabetes Technol Ther* 2012;14:411-7.

8. Zisser HC. The OmniPod Insulin Management System: the latest innovation in insulin pump therapy. *Diabetes Ther* 2010;1:10-24.
9. American Diabetes Association. 6. Glycemic targets: standards of medical care in diabetes-2021. *Diabetes Care* 2021;44(Suppl 1):S73-84.
10. Wheeler BJ, Heels K, Donaghue KC, Reith DM, Ambler GR. Insulin pump-associated adverse events in children and adolescents: a prospective study. *Diabetes Technol Ther* 2014;16:558-62.
11. Layne JE, Parkin CG, Zisser H. Efficacy of the omnipod insulin management system on glycemic control in patients with type 1 diabetes previously treated with multiple daily injections or continuous subcutaneous insulin infusion. *J Diabetes Sci Technol* 2016;10:1130-5.
12. Rice D, Sweeney K. Choosing and using an insulin pump infusion set. *Diabetes Self Manag* 2006;23:60, 2-4, 7.
13. Heinemann L. Insulin infusion sets: a critical reappraisal. *Diabetes Technol Ther* 2016;18:327-33.
14. Wang X, Cheng C, Wang S, Liu S. Electroosmotic pumps and their applications in microfluidic systems. *Microfluid Nanofluidics* 2009;6:145.
15. Beck RW, Bergenstal RM, Cheng P, Kollman C, Carlson AL, Johnson ML, et al. The relationships between time in range, hyperglycemia metrics, and HbA1c. *J Diabetes Sci Technol* 2019;13:614-26.
16. Simmons JH, Chen V, Miller KM, McGill JB, Bergenstal RM, Goland RS, et al. Differences in the management of type 1 diabetes among adults under excellent control compared with those under poor control in the T1D Exchange Clinic Registry. *Diabetes Care* 2013;36:3573-7.
17. Foster NC, Beck RW, Miller KM, Clements MA, Rickels MR, DiMeglio LA, et al. State of type 1 diabetes management and outcomes from the T1D exchange in 2016-2018. *Diabetes Technol Ther* 2019;21:66-72.
18. Jin SM. Stepwise approach to problematic hypoglycemia in Korea: educational, technological, and transplant interventions. *Endocrinol Metab (Seoul)* 2017;32:190-4.

Appendix 1. EOPatch user satisfaction survey

Research Participant's Registration Number:	Date of Survey: / / (MM-DD-YYYY)
---	----------------------------------

1. What was your previous method of injecting insulin?

- ① Multiple Daily Injection (MDI) ② Insulin pump ③ Continuous Subcutaneous Insulin Infusion (CSII)

2. Compared to conventional insulin injection methods, were you satisfied with the design, size, and weight of the ADM in this study?

- ① Very Satisfied ② Satisfied ③ Somewhat Satisfied ④ Dissatisfied ⑤ Very Dissatisfied

3. Compared to conventional insulin injection methods, do you think the menu, directions/explanations, icons, and images in the ADM were easy to understand?

- ① Very Easy To Understand ② Requires Some Thinking ③ Very Difficult To Understand

4. Compared to conventional insulin injection methods, do you think the information and data in the ADM (such as the insulin injection graph) were helpful?

- ① Very Helpful ② Helpful ③ Somewhat ④ Not Helpful ⑤ Not Helpful At All

5. Compared to conventional insulin injection methods, were you satisfied with the design, size, and weight of the patch pumps in this study?

- ① Very Satisfied ② Satisfied ③ Somewhat Satisfied ④ Dissatisfied ⑤ Very Dissatisfied

6. Compared to conventional insulin injection methods, where (on your body) did you most often apply/attach the patch pump and why?

1st: _____ 2nd: _____ 3rd: _____

Reasons: _____

7. Compared to the conventional insulin injection methods, did you find the process of attaching the patch pump to be easy/simple?

- ① Very Easy ② Easy ③ Not Easy Nor Difficult ④ Difficult ⑤ Very Difficult

8. Compared to conventional insulin injection methods, are you satisfied with your daily life while using the patch pump in this study?

- ① Very Satisfied ② Satisfied ③ Somewhat Satisfied ④ Dissatisfied ⑤ Very Dissatisfied

9. Compared to conventional insulin injection methods, what did you find inconvenient while using the patch pump for a month?

Reasons: _____

10. Compared to conventional insulin injection methods, what is your overall satisfaction level with the patch pump?

- ① Very Satisfied ② Satisfied ③ Somewhat Satisfied ④ Dissatisfied ⑤ Very Dissatisfied

11. How satisfied are you with decrease in worry about diabetes you achieved with the insulin patch pump?

- ① Very Satisfied ② Satisfied ③ Somewhat Satisfied ④ Dissatisfied ⑤ Very Dissatisfied

12. How satisfied are you with the decrease in worry about insulin injection pen/needles that you achieved with the insulin patch pump?

- ① Very Satisfied ② Satisfied ③ Somewhat Satisfied ④ Dissatisfied ⑤ Very Dissatisfied

13. What are your thoughts regarding minor problems/errors that occurred while using the insulin patch pump?

- ① Very Satisfied ② Satisfied ③ Somewhat Satisfied ④ Dissatisfied ⑤ Very Dissatisfied

14. Compared to conventional insulin injection methods, what are your opinions on the pros and cons of the patch pump in this study?

Pros: _____

Cons: _____

15. If you have any thoughts on a function that should be added to the updated product, please list them.

Reasons: _____

Research Subject's Signature and Date: _____

Research Assistant/Manager's Signature and Date: _____